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**(GB)**

**TECHNOBASE® 6000/8000**

Universal resins for veterinary medicine

**User instructions**

**Material**

TECHNOBASE 6000 and TECHNOBASE 8000 are rapid-curing resins on a methyl methacrylate basis supplied as powder and liquid.

**Indication**

Claw treatment:

Immobilisation of claws during inflammatory diseases.
Immobilization of claw tips.

Treatment of bone fractures in large and small animals:

Extra-cutaneous osteosynthesis splinting after Becker, transversal fixation by means of Kirschner drill wires, radius and ulna fractures, tibia-fibula fractures, pedal bone fractures, strengthening of dressings.

**Properties**

**TECHNOBASE 6000** is for 3 phase processing (for casting, filling and modelling).
**TECHNOBASE 8000** processes 1-phase and can be modelled directly after mixing
When it comes to modelling TECHNOBASE does not become sticky and can be modelled to any form desired without the use of aids.
The curing time is about 4 - 6 min. according to mixing ratio and product temperature.
Once cured TECHNOBASE 6000 has a non-sticky surface.
TECHNOBASE 6000 is resistant to acids, alkalis and external impurities and can be dissolved in acetone.
In a fluid and dough-y state it attacks rubber.
Polyethylene and Hostaphan films are not attacked, and can therefore be used to insulate and cover the dough.
After completion of polymerisation TECHNOBASE 6000/8000 becomes hard and abrasion-resistant allowing for mechanic machining (grinding, sawing, drilling, mil-ling etc.)

**Processing**

**Preparation**

The surface must be roughened (for example, using a hoof rasp, hoof knife or abrasive disc).

The parts getting in touch with TECHNOBASE 6000/8000 (claw, wooden block) must be free of dust and grease and completely dry (we recommend e.g. spirit or alcohol).
To achieve reasonable curing times of approx. 4 - 6 min., it may prove practical to cool the mixing vessel and material in summer or to heat it (to about 20-23°C) in winter.

**Mixing TECHNOBASE 6000 – Mixing ratio of powder/fluid about 80g:40ml (2:1)**
For batching use the beakers we have supplied.
In one beaker you batch the liquid (see graduation 40ml) and in another you batch the powder (see graduation 80g).Then pour the liquid into the powder.
Mix directly and well.The mixing ratio can be varied according to individual requirements:
Add more powder to obtain a more solid consistency and more fluid for a thinner mixture.
Directly after mixing TECHNOBASE 6000 is still fluid enough to be poured.
Then the resin turns into a dough-like consi-stency and becomes ready for filling, and a bit later, for modelling.

**Mixing TECHNOBASE 8000 – Mixing ratio of powder/fluid about 70g:25ml (3:1)**
For batching use the beakers we have supplied.
In one beaker you batch the liquid (see graduation 25ml) and in another you batch the powder (see graduation 70g).Then pour the liquid into the powder.
Mix directly and well.
When dosing individually make sure not to add too much powder.
Otherwise the paste will become too dry and cannot be modelled properly.

**Application**

Spread a thick layer of TECHNOBASE 6000/8000 on the claw sole and wood block.
Having pressed

the wood block hard onto the claw sole the surplus resin has to be modelled over the sides of the claw and the wood block.
The circumferential groove of the wood block has to be filled properly with resin to provide good stability.
After about 4-6 minutes according to mixing ratio and product temperature TECH-NOBASE 6000/8000 will cure completely and once cured will withstand high loads.

**Curing**

The curing process begins when powder and liquid are mixed.
Curing time can be shortened by heating or prolonged by cooling, but not interrupted.
Curing brings forth heat.
The resulting polymerisation temperature depends on the quantity of TECHNOBASE 6000/8000 added.
To avoid possible tissue damages, the resin can be sprinkled with cold water until curing is completed.

The synthetic resin can be cut with a hot spatula (caution - flammable!).
Openings produced in this way can be easily and permanently closed with freshly-mixed resin.
When the healing process has been completed, the plastic will be mechanically removed with the use of appropriate tools (pliers, hammer).

**Application**

**Claw treatment:**

Claws are being immobilized during inflammatory diseases by fitting a wood block to the unaffected neighbouring claw using TECHNOBASE 6000/8000 as connector.

**Bone Treatment:**

To fix the fractured ends in their original position until healing is completed, the bone ends are connected by percutaneous bone screws (Becker) or Kirschner drill wires and individually shaped bridge splints made of TECHNOBASE 6000/8000.
Advantages:
The low weight is favorable for small domestic animals and furthermore, there is the possibility of repositioning the fragments as long as the callous hasn’t hardened yet (by separating and renewed bonding over the length necessary to ensure bone healing).

**Danger advice**

Highly flammable.
Harmful by inhalation, in contact with skin and if swallowed.
Danger of cumulative effects.
Irritating to eyes, respiratory system and skin.
May cause sensitization by skin contact.

**Safety advice**

Keep container in a well-ventilated place.
Keep away from sources of ignition- No smoking.
In case of contact with eyes rinse immediately with plenty of water and seek medical advice.
Wear suitable protective clothing and gloves.
This material and its container must be disposed of as hazardous waste.

**Storage advice**

Do not store the material above 25°C .
Store in a tightly sealed container.
Do not use the material after the expiry date.
Keep out of the reach of children.

**Field of responsibility**

Our usage instructions in spoken and written form, resulting from testi carried out by us, are made to the best of our knowledge.
However, they remain without obligation, also with reference to any trade mark rights of a third party.
This does not free you from checking the products delivered by yourself that they are fit for the intended application and use.We accept no liability with regard to the application, use and processing of the product, since such use is beyond our control remaining exclusively in your field of responsibility.
Of course we guarantee for the excellent quality of our products in conformity to our General Terms and Conditions of Delivery.

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